

BONE PIN CUTTERS



JT-21.1021: Pin Cutting Pliers Multiple Action,
4-screws FRONT cutting action for soft pin up to 4.0 mm and hard pin 3.0 mm, dia
size 8 $\frac{3}{4}$ " also available in 7", 7 $\frac{1}{2}$ " 8", 8 $\frac{1}{2}$ " 9" and 9 $\frac{1}{2}$ " with
TUNGSTEN CARBIDE (TC) Tip, Gold plated handle, LIFE TIME GUARANTEE



JT-21-1016: Pin Cutting Pliers Multiple Action, 4-screws
FRONT cutting action for soft pin up to 4.0 mm and hard pin 3.0 mm, dia
size 8" also available in 7", 7 $\frac{1}{2}$ " 8 $\frac{1}{2}$ " 9" and 9 $\frac{1}{2}$ " with
TUNGSTEN CARBIDE (TC) Tip, Gold plated handle, LIFE TIME GUARANTEE



JT-21-1014: Pin Cutting Pliers Multiple Action, 4-screws
Lateral cutting action for soft pin up to 2.5 mm and hard pin 2.0 mm, dia
Size 7" also available in 7 ½" 8", 8 ½" 9" and 9 ½" with
TUNGSTEN CARBIDE (TC) Tip, Gold plated handle, LIFE TIME GUARANTEE



JT-21-1011: Pin Cutting Pliers Multiple Action, 4-screws
Front cutting action for soft pin up to 3.5 mm and hard pin 2.5 mm, dia
Size 9" also available in 7", 7 ½" 8", 8 ½" 9 1/2" and with
TUNGSTEN CARBIDE (TC) Tip, Gold plated handle, LIFE TIME GUARANTEE



JT-21-1010: JT-21-1011: Pin Cutting Pliers Multiple Action, 4-screws
Front cutting action for soft pin up to 3.5 mm and hard pin 2.5 mm, dia
Size 9 ½" also available in 7", 7 ½" 8", 8 ½" 9" and with
TUNGSTEN CARBIDE (TC) Tip, Gold plated handle, LIFE TIME GUARANTEE



JT-21-1009: Pin Cutting Pliers, Multiple Action, 4-screws,
for soft pin up to 2.00 mm and hard pin up to 1.5 mm dia.
Size 7" also available in 7 ½" 8", 8 ½" 9" and 9 ½" with
TUNGSTEN CARBIDE (TC) Tip, Gold plated handle, LIFE TIME GUARANTEE



MS-019313 Bone Pin cutters with REMOVEABLE HANDLE
Heavy duty, TC side cutting cuts pins up to 6 mm length with handle
Total length 19" with handle for easy sterilization removed 14 ½" (TC) Tip



MS-21313 Pin Cutting Pliers Multiple Action,
Heavy duty, TC side cutting OFF KIRSCHNER WIRES,
TIBIA HEAD SCREWS, BOLTS, PINS, UP TO 6 MM DIAMETER

WIRE TWISTERS



MS-1065 TC Jaw and Gold plated rings



Wire Twister 155 mm



MS-786-A
Holding forceps for cerclage wire,
size 175 mm made of stainless steel with
TC Tip and offer lifetime Guarantee



MS-786Holding forceps for cerclage wire,
Size 155 cm Made of stainless steel with
TC Tip and offer Lifetime Guarantee



M-313 Wire Twister and Cutter TC Jaws 165mm



Wire Passer 3 Sizes per set of 3-pieces

Flat nose pliers s/s



9. 1.5/2.0/2.4MM Hex Screwdriver with Sleeve (fibre Handle)
10. 2.7/3.5/4.0MM Hex Screwdriver with Sleeve (fibre Handle)



2.7mm and 3.5mm AO Bending Levers Pair 140mm



Depth Gauge



HOHMANN RETRACTOR,